

CPAP | Total Face Mask
Type 837



 **HSINER**[®]
Touching lives, Improving life.



CPAP Total Face Mask

CPAP | Total Face Mask

Type 837

- A. The transparent cover creates a crystal clear view for the user.
- B. Silicone total full face cushion fits just right and firmly on the whole face, has less pressure for the user.
- C. Anti-asphyxia valve, automatically open when the generator pressure fails, enabling patients to breathe by themselves.
- D. Quick release adaptor, press make user release / attach the tubing easily during temporary leaving.
- E. Sampling port allows measuring pressure or ETCO2.
- F. With nasogastric (NG) port.
- G. Quick release buckles, easy for user to release and take off the mask after sleep.
- H. 360° rotational elbow providing a comfortable tube position.
- I. Applicable to the users who have facial hair and general facial irregularity.

Optional Elbow Connectors



Vented
+
AAV



Non-vented



With Nasogastric(NG) Port

Total Face Mask

#13843 Small - S

#13844 Medium - M

#13845 Large - L

 **HSINER**[®]
Touching lives, Improving life.



HSINER Co., Ltd.

No. 312, Jhongshan Rd., Shengang Dist., Taichung City 429, Taiwan.

Tel: 886-4-2515-2480 Fax: 886-4-2515-2482

<http://www.hsiner.com> E-mail: sales@hsiner.com

